



# GenuSim

Project by: Ashank Behara, Conor Broeking, Yan Luo, Niviru Wijayaratne



# Problem

- Many doctors are unable to diagnose knee injuries due to a lack of exposure and experience
- Medical students need to be able to practice and learn how to diagnose knee injuries before working on human patients because practicing on injury actors is insufficient
- Lack of simulation of knee injuries in medical teaching institutions and training facilities


# Many Cases

ACL: 200,000 injuries

Meniscus Tear: >200,000 injuries

MCL: 74,000 injuries

Septic infection: >200,000 cases





## Solution


GenuSim is a knee-ligament injury simulator to help new medical students with diagnosing knee injuries





# Injuries Simulated by GenuSim

- ACL tear
- MCL tear
- PCL tear
- LCL tear
- Meniscus tear

# Product



## GenuSim vs Real Knee






## Public Link

<https://a360.co/36Nh2sT>



# CAD Model





## How GenuSim Works

- This product is a skeleton/mannequin and is not implanted in humans
- Harmonic Drive motors are positioned at the top of the femur and are connected to cables; when driven, they increase tension in cables
- **Tight = Healthy and Loose = Torn**
  - Cables simulate tension from ligaments
  - Motors drive to multiple positions to simulate a gradient of tears



# Controlling GenuSim

Each GenuSim model has a model number assigned with it. There is a **GenuSim web app** that will adjust the motors on GenuSim to simulate different injuries. By implementing this, multiple injuries can be automatically simulated at the same time with a simple click of the mouse.

# GenuSim Controller

Ashank Behara, Conor Broeking, Yan Luo, Niviru Wijayaratne

## Injury Simulation Controller

This is the web controller for the GenuSim device with model number **00000001** . On each card select the injury you want to simulate. Multiple injuries can be replicated at the same time.

### MCL Injury Adjustment

Selecting buttons will adjust the motors to tighten or loosen the strings on the simulator to mimic ligament health and injury.

Healthy

Partial Tear

Complete Tear

### ACL Injury Adjustment

Selecting buttons will adjust the motors to tighten or loosen the strings on the simulator to mimic ligament health and injury.

Healthy

Partial Tear

Complete Tear

### Septic Infection Adjustment

Selecting buttons will adjust heat at infected areas to mimic symptoms of septic infections

Healthy

Early Stage Infection

Late Stage Infection

### Meniscus Injury Adjustment

Selecting buttons will adjust the motors to tighten or loosen the strings on the simulator to mimic cartilage health and injury.

Healthy

Partial Tear

Complete Tear



## Enhancements to GenuSim

- Create an outer casing to the bone structure to mimic muscle and fat on the leg
- Pump fluids into pockets around the patella to simulate swelling when pressure is applied
- Pass heated fluids into pockets around the whole leg to simulate septic infections
- Easier accessible bones to replace with different materials i.e. replace patella with “rice crispy” patella for fractures



# Enhancements to GenuSim Web Application

- Have a feedback system which guides the student during simulation
- Create automated proficiency test
  - Changes injury every 10 minutes on the actual device
  - Enter the correct injury to get points
  - If 100% of your diagnoses are correct you are GenuSim certified
- Create authentication for users so there is one access to each model
- Create online ordering store and homepage



# Conclusion

Thank you for your time and if you have any questions we would love to answer!